

S.P.L.A.T ISA BOARD IO-953 MANUAL

Brookhaven National Laboratory
Version 1.0


BY Jack Fried (email jfried@bnl.gov)
Web page (<http://www.inst.bnl.gov/~jfried>)

INTRODUCTION

The objective behind the SPLAT ISA board was to create a timing and control system for an aerosol particle mass spectroscopy experiment. The board will be used to measure the time of flight of particles in air to trigger an ablation laser. To accurately measure the time of flight a high-speed clock with a resolution of 25ns is used. All of the boards functions are implemented in an Altera 30K PLD, which allows the board to be modified fairly simply. At present the PLD contains integrators, counters, multipliers and a PC ISA interface to implement the required function.

BOARD FEATURES

- 40 MHz Clock
- 10 MHz Clock
- Seven I/O ports
- Removable Prototype Area
- PC ISA Interface
- PC Interrupt Compatible
- High Speed PLD
- VC++ Software Interface
- Single Event Mode
- Continuous Event Mode
- Fully controllable registers



SPLAT ISA BOARD
Block Diagram

SPECIFICATIONS

INTEGRATORS (PMT1 & 2)

1. Threshold count 1 – 65,5360
2. Integration window 25ns – 6.4us
3. Saturation width 25ns – 6.375us

INTEGRATOR (CHANNELTRON)

1. Threshold count 1-65,5360 (wait for integration window to end)
2. Integration window 25ns-6.4s
3. Saturation width NONE

TRIGGER OUT DELAY (PMT 1 & 2)

1. delay trigger by -200ns - 1.6384ms

TRIGGER OUT DELAY (CHANNELTRON)

1. delay trigger by 0 - 1.6384ms

MULTIPLIER


1. multiplication constant. 0 - 1.99999

LASER TIME OUT DELAY

1. Trigger Rest Delay 1.6us - 104.8576ms

TRIGGER OUT ACCURACY

1. Time between PMT1 & PMT2 ± 50 ns
2. CHANNELTRON output trigger ± 25 ns



DIP SWITCH SETTINGS SW1

PC I/O

Address (O = open C = Closed X= NOT USED)

	{ 1 2 3 4 5 678}
0x3e0	O O O O XXX
0x3c0	O O O C XXX
0x3a0	O O O C O XXX
0x380	O O O C C XXX
0x360	O O C O O XXX
0x340	O O C O C XXX
0x320	O O C C O XXX
0x300	O O C C C XXX
0x2e0	O C O O O XXX
0x2c0	O C O O C XXX
0x2a0	O C O C O XXX
0x280	O C O C C XXX
0x260	O C C O O XXX
0x240	O C C O C XXX
0x220	O C C C O XXX
0x100	C O C C C XXX
0x1e0	C O O O O XXX
0x1c0	C O O O C XXX
0x1a0	C O O C O XXX
0x180	C O O C C XXX
0x160	C O C O O XXX
0x140	C O C O C XXX
0x120	C O C C O XXX
0x100	C O C C C XXX

PC ISA REGISTERS

Board ID

BRDID **0x00**

Function :

Return the board ID.

WRITE

NONE

READ

16 bit

Photo Multiplier 1&2 Threshold Setting

PMTTRS **0x02**

Function :

This register controls the threshold for PMT1 and PMT2.

It sets the number of pulse counts the PMTS make before a trigger is generated.

WRITE
16 bit
READ
16 bit

Photo Multiplier 1&2 Integration window Setting
PMTWND 0x04

Function:
This register sets the Integration Windows for PMT1 and PMT2
RANGE 1-256 25ns - 6.4us
1 = 25ns
2 = 50ns
.
.
255 = 6.375us
256 = 6.400us

WRITE
8 bit
READ
8 bit

Photo Multiplier 1&2 saturation Setting
PMTSAT 0x06

Function:
This register sets the saturation setting for PMT1 and PMT2
RANGE 1-256 25ns - 6.4us
1 = 25ns
2 = 50ns
.
255 = 6.375us
256 = 6.400us

WRITE
8 bit
READ
8 bit

CHANNELTRON Threshold Setting
CHNTRS 0x08

Function :
This register controls the threshold for PMT1 and PMT2.
It sets the number of pulse counts the PMTS make before a trigger is generated.

WRITE
16 bit
READ
16 bit

CHANNELTRON Integration window Setting
CHNWND 0x0a

Function:
This register sets the Integration Windows for PMT1 and PMT2
RANGE 1-256 25ns - 6.4us
1 = 25ns
2 = 50ns

.
255 = 6.375us

256 = 6.400us

WRITE 8 bit
READ 8 bit

NOT USED

CHNENC 0x0c

Function :
NONE

WRITE NONE

READ NONE

TRIGGER reset delay

LTIMEOUT 0x0e

Function :
Wait time after the trigger is generated before restarting cycle.
Range 1.6us – 104.8576ms in 1.6us steps

0 = 1.6us

1 = 3.2us

2 = 4.8us

.
WRITE 16 bit
READ 16 bit

WAIT After PMT1 Trigger Before Allowing PMT2 to Trigger

T1WAIT 0x10

Function :
Wait time after PMT1 trigger is generated before allowing PMT2.
Range 25ns – 1.638ms

1 = 25ns

2 = 50ns

2 = 75ns

.
WRITE 16 bit
READ 16 bit

TIME OUT after PMT1 Triggered But no PMT2 Trigger Arrived

T2WAIT 0x12

Function :
TIME OUT after PMT1 triggered but NO PMT2 arrives
Range 25ns – 1.638ms

1 = 25ns

2 = 50ns

2 = 75ns

WRITE
16 bit

READ
16 bit

CHANNELTRON Delay output Trigger

T3WAIT 0x14

Function : CHANNELTRON Delay output Trigger
Range 0ns – 1.638ms

0 = 0
1 = 25ns
2 = 50ns
2 = 75ns

WRITE
16 bit

READ
16 bit

NOT USED

MULTSUM 0x16

Function :
NONE

WRITE
NONE

READ
NONE

MULTIPLIER CONSTANT

ERRC x18

Function : multiplier Constant
Range 0 – 1.99999

$C1 * 4096 = 16 \text{ bit value}$

Example

T1 = 5430

C1 = 0.98

$0.98 * 4096 = 4014.08 \text{ (float)}$

reg_load -> = 4014 (16 bit)

true result = 5321.4

mult result = 5321.293

WRITE
16 bit

READ
16 bit

Photo Multiplier trigger output delay

PMTDEL 0x1a

Function : CHANNELTRON Delay output Trigger
Range -175ns – 1.6378ms

1 = -175ns delay

2 = -150ns delay

3 = -125ns delay

.

7 = -25ns delay

8 = 0ns delay

9 = 25ns delay

.

.

WRITE

16 bit

READ

16 bit

USER REGISTERS

STAR **0x1c**

Function : USER REGISTERS

Bit 0 :generate trigger 0 -> 1

Bit 1 :set 0 = single event

 :set 1 = continues trigger

Bit 2 :set 0 = PMT delay disabled

 :set 1 = PMT delay enable

Bit 4 :set 1 = reset board

 :set 0 = Normal run

Bit 14: if = 0 data from PMT

 If = 1 data from CHANNELTRON

Bit 15: if = 1 trigger fired

WRITE

14 bit

READ

16 bit

EVENT DATA OUTPUT

PTRVLT **0x1e**

Function : READ the time of flight or the # of electrons

WRITE

NONE

READ

16 bit

SOFTWARE EXAMPLE:

EXAMPLE 1

```
#include "splatisa.h"
```

```
splatisa splat; /* declare the class splatisa*/
```

```
function()
```

```

{
    /* set up splat board */
    splat.reset_board();
    splat.set_reg(PMTTRS,0x100);           //Set PMT threshold to 256
    splat.set_reg(PMTWND,0x150);         //Set window to 3.75us
    splat.set_reg(LTIMEOUT,12500)       // Set timeout to 20ms
        /*collect data */

    OPEN FILE TO SAVE DATA

    splat.cont_trigger();                 // Have the board generate continuous triggers
    while(!splat.check_trigger() && !STOP)
    {
        Save to file ← splat.read_trvl(); // Save data collected from splat
        Save to file ← splat.reg.start(); // Save what generated the trigger
        Save to file ← PCI DATA         // read data from PCI card
    }

    CLOSE FILE
}

```

EXAMPLE 2

```
#include "splat.h"
```

```
splat splat; /* declare the class splat */
```

```
function()
```

```

{
    /* set up splat board */
    splat.reset_board();
    splat.set_reg(PMTTRS,0x100);           //Set PMT threshold to 256
    splat.set_reg(PMTWND,0x150);         //Set window to 3.75us
        /*collect data */

    OPEN FILE TO SAVE DATA

    splat.single_trigger();               // Have the board generate single triggers
    while(!splat.check_trigger() && !STOP)
    {
        Save to file ← splat.read_trvl(); // Save data collected from splat
        Save to file ← splat.reg.start(); // Save what generated the trigger
        Save to file ← PCI DATA         // read data from PCI card
        splat.single_trigger();           // generate a single triggers
    }
}

```

} CLOSE FILE